

2018-08-13

Klinisk farmakologi

Cut off-nivåer

	Urin		Utandning	Saliv
	Screening	Verifikation		
Bensodiazepiner	200 ng/mL	50 ng/mL	18 pg/filter	4 ng/mL
Buprenorfin	5 ng/mL	5 ng/mL		
Diazepam	200 ng/mL	50 ng/mL	18 pg/filter	4 ng/mL
Oxazepam	200 ng/mL	50 ng/mL	18 pg/filter	4 ng/mL
Cannabis (THC-COOH)*	25 ng/mL	10 ng/mL*	18 pg/filter	4 ng/mL
Amfetamin	500 ng/mL	200 ng/mL	18 pg/filter	4 ng/mL
Metamfetamin	500 ng/mL	200 ng/mL	18 pg/filter	4 ng/mL
Ecstasy (MDMA)	500 ng/mL	200 ng/mL		
Dextropropoxifen	300 ng/mL	150 ng/mL		
Nordextropropoxifen		150 ng/mL		
Etanol	5 mmol/L	5 mmol/L		
Etylglukuronid	0,5 µg/mL	0,5 µg/mL		
Kokain (bensoylekgonin)*	150 ng/mL	100 ng/mL*	18 pg/filter	4 ng/mL
LSD	0,5 ng/mL	0,5 ng/mL		
Metadon	300 ng/mL	150 ng/mL		
Opiater	300 ng/mL			
Morfin		300 ng/mL	18 pg/filter	4 ng/mL
6-AM (6-Acetylmorfin)		10 ng/mL	18 pg/filter	4 ng/mL
Kodein		300 ng/mL	18 pg/filter	4 ng/mL
Etylmorfin		50 ng/mL		
Tramadol	200 ng/mL	50 ng/mL		
Zolpidem	50 ng/mL	50 ng/mL		
Zopiklon	50 ng/mL	50 ng/mL		